



# Shabbat



Every Amigo  
Shabbat is inspected  
and certified by a  
**Zomet Institute**  
representative.

|                            |
|----------------------------|
| <b>SHABBAT CONTROL</b>     |
| Certified by               |
| <b>The ZOMET Institute</b> |
| Approval #: _____          |
| Type: _____                |
| Model: _____               |
| Serial #: _____            |



**Halachically-approved  
for use on the Sabbath**



# Shabbat

The world's only halachically-approved mobility vehicle



- Rear drive
- Weight capacity up to 500 lbs.
- Battery capacity up to 35 miles
- 360° seat rotation
- LED battery gauge

**We go where you go...  
even on Shabbat!**

**For more information, contact:**

## How it works:

The toggle switch changes operation from "normal" to "Shabbat" mode. With no throttle activation by the rider, the module's timing circuit sets the chair in motion after approximately 10 seconds, satisfying the "no work" Sabbath restriction.

