



Seat swivels 180° for easy transfer



Headlight for increased visibility



Breaks down to 5 pieces for easy transport



Comes with rugged 9"x3" wheels

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
PHOENIXLT4	Phoenix LT 4-Wheel Travel POV 18" Seat, Grey/Blue/Red	822383012957	1/ea
PHOENIXLT4-20	Phoenix LT 4-Wheel Travel POV 20" Seat, Grey/Blue/Red	822383012964	1/ea

PRODUCT SPECS	PHOENIXLT4	PHOENIXLT4-20
Weight Capacity	350 lbs	350 lbs
Base Weight	127 lbs	130 lbs
Battery Weight (pair)	31 lbs	31 lbs
Seat Weight	23.5 lbs	26.5 lbs
Dimension (LxWxH)	42.3" (L) x 20.1" (W) x 36.6" (H)	42.3" (L) x 20.1" (W) x 36.6" (H)
Wheel Dimensions	9" x 3"	9" x 3"
Seat Dimensions	18" x 17"	20" x 16"
Top Speed	4 mph	4 mph
Maximum Range	14 miles	14 miles
Turning Radius	70"	70"
Climbing Angle	9°	9°
Ground Clearance	3.9"	3.9"
Batteries	22AH	22AH
Charger	24V/2A	24V/2A
Brakes	Electromagnetic	Electromagnetic
Motor	24V *250W x170	24V *250W x170
Controller	PG 45A	PG 45A
Warranty - Battery	12mo	12mo
Warranty - Electronics/Drive	24mo	24 mo
Warranty - Frame	Limited Lifetime*	Limited Lifetime*

BOX SPECS

Box Dimensions	49" (L) x 21" (W) x 21.5" (H)
Box Weight	143 lb

* Refer to instruction guide or website for warranty details.