



HEIGHT-ADJUSTABLE SWIVEL SEAT WITH FOLD-DOWN BACKREST

LARGE, PLASTIC CARRY BASKET

INTERCHANGEABLE BOLD COLOR PANELS

FLAT-FREE, NON-MARKING TIRES

ANTI-TIP WHEELS

EXTENDED BATTERY OPTION



Seat swivels 180° for easy transfer



Flip back armrests for easy transfer



Breaks down to 5 pieces for easy transport



Comes with blue interchangeable shroud decorations

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
SCOUTLT4	Scout LT 4-Wheel Scooter, Red/Blue	822383012933	1/ea
SCOUTLT4-EXT	Scout LT 4-Wheel 20AH Batter, Red/Blue	822383012940	1/ea

PRODUCT SPECS	SCOUTLT4	SCOUTLT4-EXT
Weight Capacity	300 lbs	300 lbs
Base Weight	100 lbs	125 lbs
Battery Weight (pair)	20 lbs	31 lbs
Seat Weight	16.6 lbs	16.6 lbs
Dimension (LxWxH)	42.3" (L) x 19.7" (W) x 36.2" (H)	42.3" (L) x 19.7" (W) x 36.2" (H)
Wheel Dimensions	8" x 2"	8" x 2"
Seat Dimensions	17.3" x 14"	17.3" x 14"
Top Speed	4 mph	4 mph
Maximum Range	9 miles	14 miles
Turning Radius	53"	53"
Climbing Angle	6°	6°
Ground Clearance	3.3"	3.3"
Batteries	12AH	20AH
Charger	24V/2A	24V/2A
Brakes	Electromagnetic	Electromagnetic
Motor	24V *250W x170	24V *250W x170
Controller	PG 45A	PG 45A
Warranty - Battery	12mo	12mo
Warranty - Electronics/Drive	24mo	24 mo
Warranty - Frame	Limited Lifetime*	Limited Lifetime*

BOX SPECS	SCOUTLT4   SCOUTLT4-EXT
Box Dimensions	46.85" (L) x 20.5" (W) x 20.5" (H)
Box Weight	114.4 lb   125.4 lb

\* Refer to instruction guide or website for warranty details.