



INVACARE® TOP END®
FORCE™ G
HANDCYCLE



The **Invacare® Top End® Force™ G Handcycle** design offers an adjustable, trunk-powered, erect position for low paras, amputees and competitors who can “gut” it out using abdominal, arm, shoulder and chest muscles. The Invacare Top End Force G Handcycle features an internal rib center frame design to make it super stiff. The adjustable back angle (55°-90°), back height, leg rests, crank height and fore/aft seat position allows you to dial in your ultimate power position. Are you ready to force the pace?

Features

Custom skeletal reinforced aluminum frame	Shimano® 27 speed drive train	Flip level parking brake
Adjustable seat length	Shimano Rapid Fire™ hands-on shifters brake	Seat and back cushions
Adjustable back angle	11-34 SRAM® cassette	Chain guard
Aluminum oversized rake fork	Cable driven upper derailleur shifter	Safety flag
Travel ready front fork	Upper chainrings: 48-39-24	250 lb. weight capacity
Adjustable crank housing	Steering dampener	
Multi-position adjustable footrests	26" high performance wheels	
Safety straps	High pressure clincher tires	
Invacare® Top End® V or S crankset	Threaded axles	
Ovalized aluminum ergonomic handpedals	Drafting bumper	

Model #: FRCG



1. Adjustable back angle 55°-90°
2. Drafting bumper
3. Cable operated crank mounted shifter for upper chainrings
4. Shimano® Rapid Fire™ Shifter/brake mounted on a handpedal with optional power plate
5. Fork with diagonal and adjustable height and chain guard
6. Travel ready front fork



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2016

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Shimano and Rapid Fire are registered trademarks of Shimano Inc. SRM is a registered trademark of SRM Service Center Inc. USA. SRAM is a registered trademark of SRAM Corp. Continental is a registered trademark of Continental Tire Corp. Specifications are subject to change without notification.
Form No. 12-450
Rev. 03/16
151173

Options

Carbon fiber wheels	Alignment gauge	Click straps
SRM® Training system	Backpack hydration system	Crutch holder
SRAM® Grip or Trigger shifters	Disk brake	Bike rack
SRAM Cassette choices: 11-23 or 11-26	Large or small hand pedals	Indoor training roller
53-40-28 Chainring	Power plates	Extra cable kit
Continental® Grand Prix clincher tires	Water bottle and cage	Leg guard attachment
Continental Sprinter tubular tires	Cordless computer	Extra crank arms
Heart rate monitor	Helmet	Tool and tire repair kit
	Safety lights, mirror	Wheel bag
	Additional safety flag	Top End bike jersey