

# ROBOOTER

R E B O O T Y O U R L I F E

## ROBOOTER E40

### A SMART WHEELCHAIR

The ultimate experience has arrived



[www.roobooter.com](http://www.roobooter.com)

# ROBOOTER E40

ULTIMATE ADVENTURE  
DELIGHTFUL INNOVATION  
YOUR ROBOOTER E40



## ENDURANCE YOU DESERVE POWERFUL & AFFORDABLE

14 Mile range with 20ah battery



## ELEGANCE MEETS FUNCTIONALITY DESIGNED FOR YOU



Ergonomic comfort

Leather stitched armrest

Advanced manufacturing,  
build to last

Durable finish

## NEXT GENERATION SOFTWARE AND HARDWARE INTEGRATION

Drive modes, automatic upgrades,  
and self-monitoring alerts; drive at ease.



CLASSIC WHITE



CLASSIC GREEN

## PRODUCT PARAMETERS

Folding	Manual folding	Front Wheel Material	Honeycomb PU solid
Maximum speed	5mph	Rear Wheel Material	Pneumatic tubeless
Maximum obstacle clearance	1.6"	Maximum capacity	330lbs.
Maximum gap clearance	4"	Battery capacity	20Ah(Optional 10ah)
Maximum climbing angle	9°	Overall dimension	(L)29"x(W)24"x(H)36"
Weight without battery	47.5lbs.	Folding shape dimension	(L)27"x(w)24"x(H)17"
Battery number	1	Seat wide	18"
Battery weight (20ah)	7.5lbs.	Maximum range	14miles (7miles with 10ah battery)
Specifications of front and rear wheels	Front: 8-in. / Rear 10-in.	Color Choices	Green, Gray, Silver, Black