

Carbon Ultralight

World's lightest rollator

Exclusively from



by  **ACRE**

With Carbon Ultralight there is no excuse to stay at home

Carbon Ultralight is the world's lightest walker at 10.6 lbs. It's designed to be the perfect companion for adventures around the world, or just around the corner. It's easy to use, easy to fold and incredibly stylish.



Compact track

Seat width: 16"
Seat height: 22"

BYA100SMR: Red
BYA100SMBK: Black
BYA100SMW: White



Regular track

Seat width: 16.5"
Seat height: 24"

BYA100MDR: Red
BYA100MDBK: Black
BYA100MDW: White



Wide track

Seat width: 18.5"
Seat height: 24"

BYA100LGR: Red
BYA100LGBK: Black
BYA100LGW: White



Weekend bag

The weekend bag offers a convenient way to carry objects and necessities for a few days away from home.

BYAWKNBG



Organizer

The organizer is designed to hold all the smaller personal belongings and function as a small purse.

BYAOBGR: Small/medium
BYAOBGL: Large



Travel bag

The travel bag is a must-have for transporting the rollator during travel.

BYATVLBG



Grocery bag

The grocery bag is designed to make your shopping easy breezy.

BYAGRBG



Backrest

The backrest offers great support for the lower back while taking a rest on the seat of the rollator.

BYASMBR: Small/medium
BYALBR: Large



Cane holder

The cane holder offers a clever way to keep the cane in place while walking with the rollator.

BYACNHR: Right
BYACNHL: Left

Carbon Ultralight features

- Only 10.6 lbs, making it the world's lightest outdoor rollator
- Ergonomic, reversed handles to secure a good posture and stability
- 8" soft and comfortable tires on all 4 wheels
- Easy height adjustment in 1" increments by just a click on a button
- Easy folding for storage and transport
- Height: min. 31", max. 36"
- Comes fully assembled

Extra warranty

By registering the product and frame number on www.byacre.com we offer an additional 3 year warranty on the frame in addition to the regular 2 years.

This warranty only applies to the frame and does not cover normal wear and moving parts.



For more information, please contact your Medline representative or visit www.medline.com/go/byacre



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093
Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996 | medline.ca | canada@medline.com
Medline México
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    